


EASY BONE SCRAPER

ZENITON 

FUNCTIONS

- The front part is bent so that the lateral access to the oral cavity is easy.
- The blade can be used 4 times in total if it is used respectively by rotating 90° as it can be rotated 360°.
- As the product is transparent, you can see in real time the amount of autologous bone that has been taken.

CODE NO. 3500



Chamber

As the chamber is made of transparent material, it can be positioned easily in the mouth and the amount of bone can be easily identified in real time.

Blade

As the blade is heat-treated, it does not easily become dull.

Body

The body is designed with the optimum angles considering ergonomics and performance.



ZENITON 